

D1008 Series

Palm Micro Centrifuge

Features

- Ideal for microfiltration and quick spin-down application
- Easy to operate. Lid close, it begin to rotate
- Max up to 7000rpm
- Compatible with 2mLx8 and PCR8x4 rotors
- Quiet running, noise level $\leq 45\text{dB}$
- Quick-change rotor without any tool
- Accelerates and brakes in seconds
- Available in two colors i.e., green, blue
- Safe, durable and robust rotor with clamp locking design.



Quiet operation



Vibration free



CE  FCC

Specifications

	D1008
Max. Speed	7000rpm
Max. RCF	2680×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×32 PCR strips or 0.2mL×4 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC100-240V/50Hz/60Hz 20W
Noise Level	$\leq 45\text{dB}$
Dimension [D×W×H]	160x170x122mm
Weight	0.5kg
Certification	CE cTUVus FCC MCA



Available rotors



A8-2P

Cat. No. 19400015


Max. speed:7000rpm

Max. RCF:2680xg

Rotor capacity:2mL/1.5mLx8,0.2mLx8,0.5mLx8

Rotor material:High strength plastic

Available rotor adapters:SA02P2, SA05P2



A04-RCR8

Cat. No. 19200810

Max. speed:7000rpm

Max. RCF:2680xg

Rotor capacity: 0.2 mlx32 PCR strips, 0.2 mlx4 PCR 8 strips

Rotor material: High strength plastic

Available adapters



SA02P2

Cat. No. 19400029

0.2mL rotor adapter

Used with A08-2 rotor

8pcs/pk



SA05P2

Cat. No. 19400026

0.5mL rotor adapter

Used with A08-2 rotor

8pcs/pk